

# NOTE TECHNIQUE

CICC

N°002 – v001

## ANTIBODY TITRATION

Why titrate antibodies ?

- 1- To **OPTIMIZE SEPARATION** of positive and negative populations
- 2- To contribute to results **CONSISTENCY AND ROBUSTNESS**
- 3- To **SAVE MONEY**.

General recommendations:

- Titrate the antibody under the conditions you will be using them.
- Set up the PMT voltages so that the negative and positive population fall in the linearity range of the detector.
- Test 2-3 concentrations around the data sheet recommendation.

« Saturating » concentration	😊	« Separating » concentration
++ for markers where expression is indicative of function or identity.		OK for markers whose expression is ON/OFF and are used as identity markers
✓ Antigen density per cell measurable ✓ Cell number should not impact results.		Level of antigen expression cannot be quantified.


1- On **HOMOGENEOUS POPULATION** (clear pos/neg peaks), calculate the stain index

$$SI = (Med_p - Med_N) / 2 * SD_N$$

Or alternatively, the Signal/ noise ratio

$$S/N = Med_p / Med_N$$

2- For populations with **HIGH VARIANCE** in fluorescence intensity, graph the MFI vs Ab concentration.



### What is the optimal concentration ?

- ✓ Maximum stain index of signal/ noise ratio.
- ✓ Point where MFI starts to plateau.
- ✓ Point where background fluorescence is comparable to unstained sample.

+ Download the Excel file for a complete tutorial and ready to use worksheets.