


The FlexMAP® 3D™ System

Staff in in charge: K. Chaba

The Next Dimension in Multiplexing



Présentation

The FlexMAP 3D system enables precise, rapid multiplexing of up to 500 unique analytes within a single sample. The dual sample fluidics paths and increased syringe injection speed facilitates a faster sample injection rate, thus allowing a very rapid acquisition of data

- Multiplex up to 500 analytes per well
- 96 and 384 well formats
- Automation compatible
- Bead-based detection of diverse analyte types, including proteins, antigens, nucleic acids
- Commercially available perconfigured kit panels

Examples of assays

BIO-RAD

Available Human Inflammation Assay Formats
Premixed all-in-one kits

Human Inflammation Panel 1, 37-Plex
(171-AL001M)

Human Inflammation Panel 1, 24-Plex
(171-AL002M)

Human Treg Cytokine Panel, 12-Plex
(171-AL003M)

